Plum Eye Wash Page 1 of 3

SAFETY DATA SHEET Plum Eye Wash

1. Identification of the substance/preparation and of the company/undertaking

Date issued	15.01.2009
Revision date	10.08.2012
Product name	Plum Eye Wash
Index no.	Version 1.5
Function	Dinging of foreign

Function Rinsing of foreign bodies and chemicals from the eyes.

Product group 4600
Use of the substance/preparation Eye wash

Manufacturer

Plum A/S
Frederik Plums Vej 2
5610
ASSENS
DANMARK
+45 6471 2112
+45 6471 2125
info@plum.dk
http://www.plum.dk
Randi Boldsen Fisker

2. Hazards identification

Description of hazard

The product is a medical device covered by Council Directive 93/42/EEC on medical devices. Consequently, the product does not require a safety data sheet and therefore this data sheet does not fully meet the statutory requirements for a safety data sheet.

3. Composition/information on ingredients			
Component name	Identification	Classification	Contents
Aqua	CAS no.: 7732-18-5		99.1 %
Sodium chloride	CAS no.: 7647-14-5 EC no.: 231-598-3		0,9 %
Column headings	European inventory of Exist name = Name as specified included in the substance lis	cts Service; EU (Einecs or Eli ing Commercial Chemical Su in the substance list (substan- st must be translated, if possil bl/vol, mg/m3, ppb, ppm, weig	bstances; Ingredient ces that are not ole). Contents given
HH/HF/HE		C = Corrosive, Xn = Harmful, F+ = Extremly flammable, F =	<u> </u>

4. First-aid measures	
General	The product is intended for first aid. If defects or adverse side effects are found during use, please report such defects or side effects to Plum A/S, info@plum.dk pr by phone +45 6471 2112, so that the necessary remedies can be initiated and the obligation to report serious adverse effects according to Appendix VII item 4 of the Directive can be observed. Serious adverse

Plum Eye Wash Page 2 of 3

effects should always be reported to the Board of Health in the state in which the side effect has occurred.

5. Fire fighting measures

Suitable extinguishing media All media can be used.

Accidental release measures

General measures No special steps are necessary.

7. Handling and storage

Handling	To ensure that the product remains sterile, the sealing must not be broken
	before use.
Storage	Best stored at 5-35°C

otorage Desiratored at 5-55 C

Special risks and properties Risk of infection from reuse of opened bottles.

8. Exposure controls/personal protection

Exposure controls

Other Information The product requires no special exposure control.

Physical and chemical properties

Physical state	Liquid
Odour	Perfume free
Colour	Transparent
Solubility description	Water-soluble
Bulk density	Value: 1002,5 kg/m ³
Boiling point / boiling range	Value: 100 °C
pH (as supplied)	Value: 5,0-7,0

10. Stability and reactivity

Stability The product is stable.

11. Toxicological information

Other information regarding health hazards

General The product has no toxicological properties.

12. Ecological information

Other ecological information

Environmental details, conclusion The product has no ecotoxicological properties.

13. Disposal considerations

Other Information There are no special circumstances relating to the disposal.

14. Transport information

Other applicable information.	The product is not dangerous goods. No special information is required for
	transport.

15. Regulatory information

julius ju	
Legislation and regulations	The product is covered by Council Directive 93/42/EEC on medical devices.

Plum Eye Wash Page 3 of 3

Consequently, the product is not covered by Regulation (EC) no 1272/2008 of the European Parliament and the Council on classification, labelling and packaging of substances and mixtures... .cf. Article 1 (5d). Regulation 1907/2006 EC of the European Parliament and Council on registration, evaluation, authorisation and restriction of chemicals (REACH) states in Article 2 (6c) that section IV regarding requirements for safety data sheets does not apply to medical equipment. Consequently, there is no requirement stating that the product mentioned above should be provided with a safety data sheet to meet the requirements of the above regulation.

16. Other information	
Further information	Expiry date for unopened, sealed bottles: At least three years from the manufacturing date. Additional information can be obtained by contacting info@plum.dk or by telephone +45 6471 2112
Information which has been added,	Changes compared to Version 1.4 (06.05.2011):
deleted or revised	Points 1, 2, 4, 7, 9, 15 og 16 are updated.
Supplier's notes	Plum Eye Wash is recommended for rinsing foreign bodies (such as dust, dirt, metal and wood splinters) and chemicals (such as solvents and oil) from the eyes. After contact with acids or bases, pH Neutral is recommended with a rinsing period of at least two minutes. Continue rinsing with Plum Eye Wash until seen by a doctor.
User notes	Turn the lid to break the seal and open the bottle. Press the bottle lightly so that the liquid will hit the eye in a soft squirt. Never press too hard.
Responsible for safety data sheet	Plum A/S